Effects of Abutment Height and Type of Cement on Bond Strength of Monolithic Zirconia Single Crowns Luted to One-Piece Zirconia Implants

Sarah Maria Blende, Dr Med Dent

Department of Prosthetic Dentistry, Center of Dentistry, Ulm University, Ulm, Germany.

Jochen Mellinghoff, Dr Med Dent Private practice, Ulm, Germany.

Sebastian Benedikt Gröller, Dr Med Dent Private practice, Rottweil, Germany.

Sigmar Schnutenhaus, Prof Dr Med Dent, MSc Private practice, Hilzingen, Germany.

Katharina Kuhn, Senior Dentist, Dr Med Dent Ralph Gunnar Luthardt, Prof Dr Med Dent Habil Heike Rudolph, Priv-Doz Dr Med Dent Department of Prosthetic Dentistry, Center of Dentistry, Ulm University, Ulm, Germany.

Purpose: To test the bond strength of one-piece zirconia implants with either standard or reduced height using different luting agents and pretreatments of the inner surfaces of ceramic crowns. *Materials and* Methods: A total of 20 monolithic zirconia CAD/CAM-fabricated single crowns per luting agent group were cemented onto 20 one-piece zirconia implants with either 5-mm (n = 10) or 4-mm (n = 10) abutment height using 13 different luting agents. After water storage, the crowns were removed using a specially developed test fixture in a universal testing machine. The maximum force (N) was recorded, and the force per area (MPa) was calculated. The statistical evaluation was performed using univariate analysis of variance. Results: A mean of 4.19 MPa (SD 2.90) at 5-mm height and 3.89 MPa (SD 2.85) at 4-mm height was obtained for all luting agents. The highest values were achieved for a resin-modified glass-ionomer cement, with 12.37 MPa (4 mm) and 12.00 MPa (5 mm). The lowest values were shown for a long-term temporary material, with 0.73 MPa (4 mm) and 1.07 MPa (5 mm). Only a polycarboxylate cement (P < .001) and a glassionomer cement (P = .006) showed statistically significant differences in favor of the reduced abutment height. Reduced abutment height did not significantly reduce bond strength for any of the materials examined. Conclusion: Implants with a reduced abutment height are clinically suitable. Pretreatment of the crowns' inner surfaces with ceramic primer showed to be advantageous. Int J Prosthodont 2022;35:XX–XX. doi: 10.11607/ijp.7110

ccording to current standards in clinical treatment, implant rehabilitation is the method of choice for replacing missing teeth unless it is contraindicated by the clinical situation and/or the patient's medical history. With growing demand for implant-based solutions in dentistry, the continuous development of new implant constructions and materials has become an important component of dental research.

Over the past decades, implants made of titanium have been considered the gold standard in the selection of materials. However, despite all their positive characteristics—like mechanical strength, biocompatibility, and clinically excellent survival rates—the gray color of the material can be challenging in esthetically relevant areas,^{1,2} as the grayish shimmering of the implant leads to dissatisfaction among patients and dentists, particularly in combination with thin overlying mucosal tissue.^{3,4} The Correspondence to: Dr Sarah Maria Blender Department of Prosthetic Dentistry Center of Dentistry, Ulm University Albert-Einstein-Allee 11, 8908 Ulm, Germany Fax: +49 731 500 64247 Email: sarah.blender@ uniklinik-ulm.de

> Submitted April 28, 2020; accepted February 16, 2021. ©2022 by Quintessence Publishing Co Inc.

accumulation of titanium particles in surrounding tissue and even in regional lymph nodes after implant placement has also been described. A correlation between the amount of particles and the time since implantation has been observed.^{5,6} Although little clinical evidence has been reported, hypersensitivity and local symptoms due to titanium implants have been described.^{1,7} Patient demand for metal-free alternatives to reduce the foreign body sensation and to optimize the esthetic outcome has also increased in recent years.³ With the development of modern high-performance ceramics, new alternative implant materials are now available. In particular, polycrystalline yttrium-stabilized zirconium oxide (3Y-TZP), with its excellent mechanical properties, is the most frequently used material for ceramic implants.⁸ In addition to the well-described biocompatibility of this material, its low affinity for plague formation plays an important role in the reduction of peri-implant inflammatory processes, thus promoting osseointegration of the implant.⁹ Most zirconia implants currently available are one-piece systems made of Y-TZP. Even though long-term study results are rare, initial clinical results are satisfactory, showing survival rates of 95.6% after 1 year¹⁰ and 95% after 5 years,¹¹ which are comparable to titanium implants. The same applies to marginal bone loss, which can also be a parameter for successful implant treatment. With an average rate of 0.79 mm of bone loss after 1 year, there is no significant difference compared to titanium implants.¹⁰ Two-piece zirconia implant systems are also available but show higher average loss rates due to a lack of osseointegration. The main reason is likely the gap between abutment and implant as a potential source of risk for inflammation.¹²

Implant-supported restorations can be either screw retained or cemented, and the retention modalities show no significant differences in survival rates after 5 years. Cemented restorations in particular have several advantages, such as a cost-effective and simple workflow, excluded loss of the nonexistent implant screw, and also improved and simplified results in esthetics and occlusion due to the absence of a screw channel.¹³ However, cemented restorations are more prone to biologic complications, as the presence of an adhesive gap increases the risk of peri-implant inflammation of the surrounding tissue. Special attention must be paid to the remaining luting material in the peri-implant sulcus after cementation, as this can lead to bone loss.^{14,15} Complete healing of the peri-implant tissue can be expected with subsequent removal of excess material.

When selecting a suitable restorative material, allceramic systems made of lithium disilicate or zirconium oxide are the materials of choice. Due to advances and material innovations in the field of CAD/CAM technology, monolithic fabricated restorations are preferred over veneered all-ceramic superstructures.¹² Monolithic restorations are usually cemented or adhesively seated due to a lack of available ceramic implant systems that allow screw fixation of restorations.

Furthermore, dentists are often challenged by a lack of vertical space, making the prosthetic restoration of natural teeth or implant abutments complex and difficult. The influence of various preparation parameters of natural teeth on the retention and survival of dental restorations has been well described in the literature.^{16,17} Precise guidelines and regulations for the height of abutments are so far missing in standard German textbooks, but the necessary height of prepared teeth is well covered. Both Reitemeier et al and Strub et al require an axial wall height of at least 3 mm after preparation in order to ensure sufficient mechanical retention of the future dental restoration.^{18,19} When reducing the axial wall height, the surface area available for bonding is minimized, which is also important for the retention of dental restorations.¹⁷ Additional accessories, such as boxes and grooves, can compensate for unfavorable preparation parameters, but are contraindicated when using CAD/CAM-fabricated restorations, as these demand strict avoidance of undercuts and sharp edges.²⁰ Therefore, the problem of adequate retention requires proper selection of the luting agent and good knowledge of the available materials and their contents to compensate for eventual deficiencies of the preparation. While bonding with conventional materials is based on mechanical friction between the surfaces and the cement, adhesive luting agents in combination with a ceramic primer can adhere directly to both surfaces (tooth/abutment and inner surface of the restoration) to increase retention.²¹ In addition to established materials, new developments should also be considered within the scope of an investigation.

The aim of this study was to evaluate the effects of reducing the abutment height of one-piece zirconia implants and to assess the influence of different luting agents and different surface-conditioning methods on the retention of monolithic zirconia single crowns.

MATERIALS AND METHODS

The specimens were divided into 16 test groups (n = 20) depending on the selected luting agent and the conditioning method for the inner surface of the ceramic crown (Table 1). In each test group, 10 zirconia one-piece implants with a standard abutment height of 5 mm and 10 zirconia one-piece implants with a reduced abutment of 4 mm were used, resulting in 320 test specimens. All implants (Zirkolith Z5m implants, Z-Systems) were made of TZP-A BIO-HIP ceramic (Metoxit). Compared to the abutments with reduced height, the standard abutments had a reduced upper diameter but an enlarged bonding surface due to

Material	Manufacturer	Type of application	Light-curing/application of glycerine gel	Pretreatment of ceramic crowns
HIS	Harvard Dental	Automix syringe	Light-curing 20 s per side	Al ₂ O ₃ 50 µm, 10 s
CIT	International		Light-curring 20 s per side	Cleaning with ethanol
DUR	ЗМ	Automix capsule		Al ₂ O ₃ 50μm, 10 s
DOIN	IVIC		_	Cleaning with ethanol
PHZ	Kulzer	Manual mixture of powder and liquid		Al ₂ O ₃ 50 μm, 10 s
1112	Kuizei			Cleaning with ethanol
		Automix syringe	Liquid Strip Transparent, Ivoclar	Al ₂ O ₃ 50 μm, 10 s
BQM	Voco		Vivadent	Ceramic Bond, Voco
			Light-curing 20 s per side	Cleaning with ethanol
DOD			Liquid Strip Transparent	Al ₂ O ₃ 50 μm, 10 s
BOP	Voco	Automix syringe	Light-curing 20 s per side	Cleaning with ethanol
121	Kungan Newitela	Manual mixture of two pastes		Al ₂ O ₃ 50 μm, 10 s
P21	Kuraray Noritake		Oxyguard II, Kuraray Noritake	Cleaning with ethanol
				Al ₂ O ₃ 50 μm, 10 s
PV5	Kuraray Noritake	Automix syringe	Light-curing 20 s per side	Clearfil Ceramic Primer Plus
				Cleaning with ethanol
PUO	Kuraray Noritake	Automix syringe	Light-curing 20 s per side	Cleaning with ethanol
PCO	Kuraray Noritako	Automix arringa	Light curing 20 c par side	Al ₂ O ₃ 50 µm, 10 s
PGO	Kuraray Noritake	Automix syringe	Light-curing 20 s per side	Cleaning with ethanol
BSE	Masa	Automic arriage	Light-curing 20 s per side	Al ₂ O ₃ 50 µm, 10 s
DSE	Voco	Automix syringe	Liquid Strip Transparent	Cleaning with ethanol
REX	ЗM	Automix capsule		Al ₂ O ₃ 50 µm, 10 s
REA.	5101	Automix capsule	_	Cleaning with ethanol
		Automix syringe	Liquid Strip Transparent	Al ₂ O ₃ 50 µm, 10 s
SCP	Ivoclar Vivadent			lvoclean, lvoclar Vivadent
				Cleaning with ethanol
PC 4		Automix syringe	-	Al ₂ O ₃ 50 μm, 10 s
PSA	Kuraray Noritake			Cleaning with ethanol
AQM		Manual mixture of powder and liquid		Al ₂ O ₃ 50 μm, 10 s
	Voco		-	Cleaning with ethanol
KEC.	ЗM	Automix capsule		Al ₂ O ₃ 50 μm, 10 s
KEC			-	Cleaning with ethanol
	GC Europe	Automix syringe		Al ₂ O ₃ 50 μm, 10 s
FUJ			-	Cleaning with ethanol

Table 1	Materials	Used and the Diff	erent Conditionin	a Methods for the	e Inner Surface of the	Ceramic Crowns

HIS = Harvard Implant Semi-Permanent; DUR = Durelon; PHZ = Phosphatzement; BQM = Bifix QM; BOP = Bifix Qm, sandblasting, no primer; P21 = Panava 21; PV5 = Panavia V5; PUO = Panavia V5, no sandblasting, no primer; PGO = Panavia V5, sandblasting, no primer; BSE = Bifix SE; REX = RelyX Unicem; SCP = SpeedCEM Plus; PSA = Panavia SA Cement Plus; AQM = Aqua Meron; KEC = Ketac Ceram; FUJ = Fujicem; .



Fig 1 Ceramic crown and one-piece zirconia implant. The ceramic crown was deliberately overcontoured (1.5–2 mm) in the marginal area to prevent the assembled implant-crown complex from sliding out of the test fixture during the debonding process in the universal testing machine.



Fig 2 Specimen fixed to the lower mounting box of the universal testing machine.

the different heights of the abutments with the same taper (7 degrees) and the same lower diameter (4.4 mm). To prevent rotation, the abutment design had two parallel surfaces with grooves on opposite sides (Fig 1). A total of 20 monolithic CAD/CAM–fabricated zirconia single-tooth crowns (Ceramill Zolid HT+ Preshades, Amann Girrbach) were produced by a dental laboratory (Joachim Maier Dental Design). Designed as a maxillary premolar, the marginal area of the crown was deliberately overcontoured (1.5 to 2 mm) to prevent the assembled implant-crown complex from sliding out of the test fixture during the debonding process in the universal testing machine (Z010, Zwick Roell; Fig 1).

Prior to cementation, all abutments were cleaned with alcohol and then dried with oil-free air. Any further pretreatment of the surface was not required by the manufacturer. The cementation area of the onepiece implants (abutment part) showed no optical or tactile surface roughness. The abutment surface was not sandblasted at any time (this applies to pretreatment for bonding as well as all necessary steps of the specimens' preparation for reuse). Before cementation, the abutment surface was cleaned with alcohol. After debonding, the cement residues were carefully removed with a laboratory cloth. In some rare cases, small remains were carefully sheared off with a scalpel, which did not damage the integrity of the abutment surface (light microscopic inspection).

The inside of the ceramic crowns was sandblasted with aluminum oxide particles (50 µm, 2 bar, 10 seconds, held at a distance of 10 mm), cleaned with ethanol in an ultrasonic bath (5 minutes), and dried with oil-free air. Depending on the luting material used, further pretreatment of the surface was applied where required (Table 1). The luting agents were mixed according to the manufacturers' instructions. Light curing and a glycerin gel were applied when indicated. The ceramic crowns were filled with the luting agent and then placed on the abutments using finger pressure, and excess marginal luting material was removed with a microbrush. All specimens were loaded with a force of 5 kg for the indicated setting time at room temperature (24°C) and then stored in a warming cabinet at 37°C in distilled water. After 24 hours, the specimens were fixed in the universal testing machine to detach the restorations from the implant abutments. For this purpose, an experimental test fixture was developed based on an established testing method for investigating the shear and tensile bond strength of dental materials,^{22,23} which are simultaneously present with the use of anatomical crowns. In order to ensure pull-off forces in the axial direction, the upper part of each specimen was fixed in the testing machine using a mounting box with a bottom opening of 4.8 mm, analogous to the thickest part of the implant (Figs 2 and 3). The overcontoured margins of the crowns prevented



that can be used to attach the mounting box with a hole in the bottom. (b) The implant part of the test specimen has a smaller diameter than the hole in the mounting box, which allows the implant part to slide through the bottom. The overcontoured crown has a larger diameter than the hole, making it impossible for the test specimen to slip through. The mounting box including the crownimplant test specimen can be connected to the upper part of the testing machine via the internal thread. (c) Only the implant part extends out of the box. (d and e) The implant part of the test specimen is then screwed into the internal thread of the lower mounting box. (f) This fixture is then connected to the universal testing machine via another two shekels.

the specimen from slipping through the hole in the bottom of the mounting box during crown removal (Fig 3b). The external threads of the implants were used to attach the implant to the lower holding device (Figs 3d to 3f).

Prior to debonding, the parameters of the universal testing machine were checked for analysis: crosshead speed of 0.5 mm/minute, preload of 5 N, a force limit of 10,000 N, and a force threshold of 20 N. The maximum

Fundamental Research

Test series	n	Bond strength	Retention force
HIS	20	0.92 ± 0.36	53.50 ± 23.08
DUR	19	1.62 ± 1.02	89.96 ± 51.20
PHZ	20	2.98 ± 1.25	170.49 ± 68.02
BQM	20	1.67 ± 0.32	95.71 ± 16.47
P21	20	3.60 ± 0.51	207.75 ± 36.51
PV5	20	4.41 ± 0.97	253.70 ± 58,84
BSE	20	2.67 ± 0.84	153,99 ± 53.17
REX	20	3.71 ± 0.92	213.40 ± 55.45
SCP	20	4.30 ± 0.72	248.80 ± 51.36
PSA	20	5.40 ± 0.85	309.60 ± 46.81
AQM	20	3.68 ± 1.22	212.10 ± 70.60
KEC	20	5,20 ± 1.13	296.95 ± 54.98
FUJ	20	12.18 ± 1.82	700.95 ± 117.23

Table 2	Mean ± SD Bond Strength	(MPa) and Retention Force (N) of the Different Luting Agents
---------	-------------------------	------------------------------	----------------------------------

See Table 1 for abbreviations.

retention force at the time of debonding was recorded in Newtons (testXpert II, Zwick Roell). To receive the combined shear bond and tensile strength in megapascals, the surface area of the abutments was calculated based on the manufacturers' construction drawings. Based on these results (mm²), the recorded maximum force was used to calculate the force per area (N/mm² = MPa).

After debonding, the crowns and implants were meticulously cleaned. First, careful mechanical removal of larger cement residue was done with a scalpel. The abutments and crowns were then checked under the light microscope for any further cement residue. In case of remaining residue, the inside of the crown was completely cleaned by sandblasting with alumina particles (25 µm, 2 bar, held at a distance of 10 mm). For reuse, the surface of the abutments was treated with alcohol, laboratory cloths, and, in a few cases, with a scalpel, since there were only a few abutment surfaces with cement residue. A light microscope (see below) was used to ensure an entirely residue-free surface.

To qualitatively assess the fracture behavior of the luting agents, a light microscopic examination (magnification ×2.5) was carried out (Leica MS5). The fracture behavior was classified into two different categories: (1) adhesive fracture between two different surfaces, therefore either between the luting agent and the abutment or between the luting agent and the inner surface of the ceramic crowns; or (2) cohesive fracture within the ceramic or luting agent. In order to obtain a quantitative evaluation, the bonding surface was divided into different areas. The percentage of the total surface and the surface in the individual areas was calculated, resulting in an exact description of the individual fracture pattern in percent. This analysis was always performed first.

Afterwards, the specimens were cleaned for reuse as described above.

The data acquired were statistically analyzed with a suitable statistical software (SPSS version 25, IBM). Subsequent to the descriptive statistics and *t* test for two independent samples to evaluate the influence of the different abutment heights on the bond strength in the different test series, analysis of variance (ANOVA) was performed to determine the influence of surface pretreatment. In order to identify which subgroups showed a significantly different outcome, Tukey post hoc test was applied. The level of significance was set at $\alpha = .05$ for all analyses.

RESULTS

All 320 implant-crown complexes were detached from each other. One specimen was excluded from data analysis due to incorrect fixation in the universal testing machine.

Influence of Different Luting Agents

The mean bond strength values (\pm SD) and the maximum retention force at the time of debonding (\pm SD) are shown in Table 2 for the 13 different luting agents independent of the different abutment heights. Materials from the group of glass-ionomer cements and self-adhesive luting agents showed the highest bond strength and retentive forces, followed by adhesive luting agents and conventional materials like polycarboxy-late cement (Durelon [DUR], 3M) and phosphate cement (Phosphatzement [PHZ], Kulzer). The lowest values were shown for temporary luting agents.

Test series	Abutment height (mm)	n	Bond strength	P value	Retention force
HIS	4	10	0.76 ± 0.27	.051	40.48 ± 14.55
	5	10	1.07 ± 0.37		66.51 ± 23.16
	4	9	2.51 ± 0.77	004	133.06 ± 41.10
DUR	5	10	0.82 ± 0.23	< .001	51.18 ± 14.55
	4	10	3.19 ± 1.59	460	168.77 ± 84.09
PHZ	5	10	2.77 ± 0.83	.468	172.20 ± 51.86
DOM	4	10	1.80 ± 0.32	067	95.40 ± 17.06
BQM	5	10	1.54 ± 0.27	.067	96.09 ± 16.77
0.2.1	4	10	3.52 ± 0.49	100	186.40 ± 25.80
P21	5	10	3.68 ± 0.54	.496	229.10 ± 33.70
	4	10	4.43 ± 1.06	.919	234.50 ± 56.31
PV5	5	10	4.38 ± 0.93		272.90 ± 57.61
DCE	4	10	2.58 ± 0.66	.644	136.37 ± 34.85
BSE	5	10	2.76 ± 1.02		171.60 ± 63.75
	4	10	3.78 ± 0.85	750	199.90 ± 44.95
REX	5	10	3.65 ± 1.03	.759	226.90 ± 63.77
SCP	4	10	4.14 ± 0.60	210	219.10 ± 31.90
	5	10	4.47 ± 0.82	.310	278.50 ± 50.90
	4	10	5.70 ± 0.82	100	301.70 ± 43.46
PSA	5	10	5.10 ± 0.82	.122	317.50 ± 50.97
0.014	4	10	3.63 ± 1.51	.856	192.00 ± 79.63
AQM	5	10	3.73 ± 0.92		232.20 ± 57.31
KEC	4	10	5.86 ± 0.85	225	310.20 ± 44.84
	5	10	4.56 ± 1.01	.006	283.70 ± 63.10
5111	4	10	12.37 ± 1.49	.658	655.00 ± 78.93
FUJ	5	10	12.00 ± 2.16		746.90 ± 134.50

Table 3Mean ± SD Bond Strength (MPa) and Retention Force (N) of the Different Luting Agents Depending
on the Different Abutment Heights

See Table 1 for abbreviations.

Influence of Different Abutment Heights

When comparing the bond strength values of all specimens with a reduced abutment height (n = 129) to abutments with a standard height (n = 130) regardless of the different luting agents, there were no significant differences in bond strength (4 mm: 4.19 ± 2.90 MPa; 5 mm: 3.89 ± 2.85 MPa) or maximum retention force at the time of debonding (4 mm: 221.67 ± 153.34 N; 5 mm: 241.94 ± 177.33 N).

The mean (\pm SD) bond strength in MPa for each test series and its subgroups are shown in Table 3. The use of a polycarboxylate cement (DUR; P < .001) and a glass-ionomer cement (Ketac Cem [KEC], 3M; P = .006) showed statistically significant differences, while all other materials showed no significant differences in bond strength when evaluating the influence of abutment height in relation to the different luting


agents. However, the statistically significant difference was in favor of the lower abutment height.


Influence of Surface Pretreatment

When using primers with adhesive resin-based luting agents (Bifix QM [BQM], Voco; without primer: 4 mm vs 5 mm: P = .917; with primer: 4 mm vs 5 mm: P = .169; and Panavia V5 [PV5], Kuraray Noritake; 4 mm [no primer] vs 4 mm [primer used]: P < .001; 5 mm [no primer] vs 5 mm [primer used]: P < .001), there was a significant difference in bond strength for both abutment heights (Fig 5). The influence of pure sandblasting when no ceramic primer was used showed no significant difference when using the adhesive resin-based luting agents PV5 (4 mm [no sandblasting or use of primer] vs 4 mm [sandblasting and no primer used]: P = .922; 5 mm Fundamental Research

Abutment height (mm) 4 mm 5 mm 1.50 0 4 mm 5 mm 5 mm BOP (without primer) BQM

Fig 4 Bond strength values (MPa) of the adhesive cement Bifix QM (BQM) using different conditioning methods for pretreatment of the inner surface of the ceramic crown.





[no sandblasting or use of primer] vs 5 mm [sandblasting and no primer used]: P = .539) (Fig 5).

Fracture Behavior Analysis

The resin-based glass-ionomer cement (FujiCem 2 [FUJ], GC Europe) specimens mainly showed cohesive fracture behavior within the luting material and small amounts of adhesive fractures. In the test series with a polycarboxylate cement (DUR), a phosphate cement (PHZ), two adhesive resin-based luting agents (BQM and Panavia 21 [P21], Kuraray Noritake), and a self-adhesive resinbased luting agent (Bifix SE [BSE], Voco), most of the luting material was found on the inner surface of the ceramic crowns, with small parts also detected on the lower half of the abutments. Examining the use of a temporary cement (Harvard Implant Semi-Permanent [HIS], Harvard Dental International), an adhesive resinbased luting agent (PV5), three adhesive resin-based luting agents (RelyX Unicem [REX], 3M; SpeedCEM Plus [SCP], Ivoclar Vivadent; and Panavia SA Cement Plus [PSA], Kuraray Noritake), and two glass-ionomer cements (Agua Meron [AQM], Voco; and KEC), most of the luting material was found on the inner surface of the crowns and only a small portion on the lower half of the abutments. Cohesive fractures could only be observed sporadically.

DISCUSSION

For the one-piece ceramic implants and the existing abutment geometry evaluated herein, no negative influence on bond strength was found when examining a reduced abutment height compared to the standard height, irrespective of the luting material or pretreatment methods.

Comparison of the study results to those of similar studies is very difficult due to the multitude of different experimental setups and parameters, and only a few studies are comparable. Even the many research reports with a similar experimental setup and similarly selected parameters examine a wide range of different materials for single-tooth crowns and abutments. Only studies with zirconia ceramic abutments and monolithic zirconia ceramic crowns evaluated with a crown pull-off test can be considered comparable to the study at hand. As long as there are no standardized test procedures for determining retention and resistance behavior using standardized test setups, even small changes in the parameters during implementation can have a major impact on the results obtained.^{24–27} Likely the biggest difference when compared to other test set-ups is the fixation of the test specimen in the universal testing machine. The specimen is often attached to the upper jig of the machine by wires or special crowns with built-in loops, while the implant or abutment is polymerized in resin or similar materials and is in contact with the lower part of the machine via an indirect connection.^{24,28–37} The abutments or implants in turn are fixed with resins or clamps to the lower part.

In the experimental setup presented here, a specially developed test fixture was used. It is a further advancement of an established fixture for shear and tensile bond strength testing.^{38,39} The type of upper and lower fixation was chosen such that the tensile load is applied purely axially without any influence of a load moment or torque. This also ensures that only the contact areas between the materials of interest (ie, the implant-crown complexes themselves) are tested while avoiding confounders due to additional fixation materials and the kind of fixation, often indirect. The upper fixation consists of a mounting box, which has a bottom opening analogous to the width of the implant and which can be fixed by an internal thread to the external thread of the upper part of the machine. The holding box is then put over the external thread of the implant but stops in the area of the crown margin because the crown walls have been thickened. In the lower area of the testing machine, the implant is fixed directly over the thread in the testing machine.

In the present study, the bond strength values of the individual luting agents showed a wide range, but a trend between the individual material classes could be observed. Comparing the achieved bond strengths of the different classes of luting materials achieves the following ranking: glass-ionomer cements > self-adhesive resin-based luting agents > resin-based adhesive luting agents > phosphate cement > polycarboxylate cement > temporary cements.

This can be attributed to the different mechanisms of attachment of the luting agents. Conventional materials are held in place by mechanical mechanisms such as friction, which is why the design of the preparation and the surface quality play an important role here.²¹ These materials therefore necessitate a very close fit between the restoration and the natural tooth/implant abutment in order to achieve sufficient retention.

In adhesive materials, a chemical component is used in addition to the macromechanical retention, which is clearly noticeable in the interaction of luting material and tooth or abutment surface and the resulting retention strength. Several studies particularly recommend the use of dual or self-adhesive resin-based luting agents when zirconia ceramics are used.^{40–43} Since removability is an inherent property of (semi)permanent luting materials, it is not surprising that the lowest adhesion values are achieved with these.

Some studies showed similar values of bond strength for the same materials tested in the present study.^{33–37} In particular, differences were found for phosphate cements when compared to other studies. These cements showed much higher bond strength values compared to the present study.³⁵ The retention force of phosphate cement is the result of friction only. As no surface treatment was applied to the implant abutments used in this study, resulting in a smooth surface and a lower retention, this smoothness could be a reason why luting agents with adhesive molecules showed to some degree higher bond strength values than conventional materials. The resin-based glass-ionomer cement showed the highest values for bond strength. Due to the strength of the adhesive bond, a direct chemical bond between the luting material and the zirconia surfaces cannot be ruled out. In the authors' opinion, this bond cannot be explained by friction alone.

The reduction of the abutment height showed significant differences in the test series when using a polycarboxylate cement (DUR) and a glass-ionomer cement (KEC).

The fact that 4-mm abutments (ie, those with reduced height) also showed significantly higher values depending on the respective luting material allows the conclusion that, although statistically significant, there is no clinical relevance of the abutment height of the onepiece implants investigated here. Otherwise, the 5-mm abutments should have consistently shown the higher bond strength values. This conclusion is based solely on the comparison of the abutment geometries tested here.

A similar in vitro study examined the influence of 4-mm and 5.5-mm zirconia abutments and found significant differences in the use of temporary cements, self-adhesive materials, and glass-ionomer cements. The possible use of an additional ceramic primer and the pretreatment of the zirconia abutments was not described.³⁵ Other studies using different materials for abutments and crowns showed that the influence of several abutments with a height of 5 mm or less did not significantly differ, while abutments with a height of 5 mm or more showed an increase in retention.

Different pretreatment methods for zirconia crowns are controversially discussed in the literature. The results of this study showed that the use of sandblasting with aluminum oxide particles and the use of ceramic primers led to an increase in bond strength. These observations are also confirmed by other studies.^{30,44–46} Roughening of the surface results in expansion and improved wettability of the contact area. When using resin-based luting agents, an additional treatment of the surface with a silane coupling agent or ceramic primer with 10-methacryloyloxydecyl dihydrogen phosphate (MDP) molecules showed significantly increased bond strength values.^{44,46}

The fracture analysis showed mainly adhesive fracture behavior and thus revealed a weak point in the bond between the respective luting agents and the smooth abutment surface. Only the luting agents with the highest bonding strength values showed areas of cohesive fracture. According to the manufacturers' instructions, additional surface pretreatment was not required. In this study, a clear improvement in bond strength could be achieved after pretreating the abutment surface with appropriate means, such as alumina particle sandblasting. Several studies showed that a significant improvement in bond strength could be reached on both natural teeth and titanium abutments when using sandblasting products for surface pretreatment of the inner surfaces of restorations.^{29,36,47}

Currently, there is no laboratory test procedure with which the clinical parameters of the oral cavity can be exactly reproduced. This should be considered when looking at in vitro results.^{24,29,48} However, laboratory tests have the advantages of being cost effective, clearly reproducible, and a rapid alternative to clinical studies, which is important for estimating the technical properties of new materials.^{24,26,29} As temperature may influence the setting or curing of luting materials, the implants were stored at 37°C until immediately before luting and removed from the warming cabinet only for the actual luting process. In addition, the setting time at body temperature specified by the manufacturer was doubled, and the connected specimens were immediately returned to the warming cabinet to take this aspect into account.

In order to be able to reuse the single-tooth crowns, cement residues were removed and the surfaces cleaned as described above. The cementation gap was not measured before the bonding of the individual test series nor after each cleaning, as other in vitro investigations of the research group using densely sintered ZrO_2 specimens showed an eventual loss of surface material in the single-digit micrometer range (less than 6 µm) after sandblasting with aluminum oxide particles of 25 µm in size and below 10-bar pressure and repeated noncontact optical measurements.⁴⁹ While a more precise measurement was not possible due to the lack of a suitable measuring method at the research site, the potential loss was assessed as negligible.

The negative influence of long-term water storage or thermocycling, the existence of which is known,^{50,51} was not the subject of this in vitro study and should be the subject of future investigations.

The present results apply only to abutments with the same outer shape and taper. In addition to abutment

height, even minimal changes in parameters such as abutment shape or taper can lead to major differences in study results. Therefore, only results of studies with similar testing procedures and parameters are comparable.^{17,24–27}

CONCLUSIONS

The proper selection of a luting agent when fixing zirconia single crowns onto one-piece zirconia implants is important. Materials classified in the group of glass-ionomer–based cements or self-adhesive resins showed the highest bond strength values, followed by adhesive resins, conventional luting agents, and temporary materials.

Depending on the use of different luting agents, only a minority of the selected materials showed significant differences in bond strength in favor of the lower abutment height, thus being inconsequential for clinical application. For this kind of one-piece zirconia implant, comparable bond strength values can be achieved with both abutment heights under in vitro conditions.

Surface pretreatment with ceramic primers showed significantly higher bond strength values.

ACKNOWLEDGMENTS

This study was supported in part by Z-Systems and Joachim Maier Dental Design. The results have been presented in part (oral presentation) at: the 16th Voco Dental Challenge, September 21, 2018, Cuxhaven, Germany (third place); the 68th Annual Meeting of the German Society for Prosthetic Dentistry and Biomaterials (DGPro), May 17–18, 2019, Rostock, Germany. The authors report no conflicts of interest.

REFERENCES

- Osman RB, Swain MV. A critical review of dental implant materials with an emphasis on titanium versus zirconia. Materials (Basel) 2015;8:932–958.
- Smith DC. Dental implants: Materials and design considerations. Int J Prosthodont 1993;6:106–117.
- Andreiotelli M, Wenz HJ, Kohal RJ. Are ceramic implants a viable alternative to titanium implants? A systematic literature review. Clin Oral Implants Res 2009;20(suppl 4):32–47.
- van Brakel R, Noordmans HJ, Frenken J, de Roode R, de Wit GC, Cune MS. The effect of zirconia and titanium implant abutments on light reflection of the supporting soft tissues. Clin Oral Implants Res 2011;22:1172–1178.
- Bianco PD, Ducheyne P, Cuckler JM. Local accumulation of titanium released from a titanium implant in the absence of wear. J Biomed Mater Res 1996;31:227–234.
- Weingart D, Steinemann S, Schilli W, et al. Titanium deposition in regional lymph nodes after insertion of titanium screw implants in maxillofacial region. Int J Oral Maxillofac Surg 1994;23:450–452.
- de Graaf NPJ, Feilzer AJ, Kleverlaan CJ, Bontkes H, Gibbs S, Rustemeyer T. A retrospective study on titanium sensitivity: Patch test materials and manifestations. Contact Dermatitis 2018;79:85–90.
- Dahane TM, Pakhan AJ, Godbole D, Kale D, Bhola R. A review on zirconia as the upcoming material for dental implants. IOSR-JDMS 2017;16:78–85.

- Hashim D, Cionca N, Courvoisier DS, Mombelli A. A systematic review of the clinical survival of zirconia implants. Clin Oral Investig 2016;20:1403–1417.
- Pieralli S, Kohal RJ, Jung RE, Vach K, Spies BC. Clinical outcomes of zirconia dental implants: A systematic review. J Dent Res 2017;96:38–46.
- Oliva J, Oliva X, Oliva JD. Five-year success rate of 831 consecutively placed Zirconia dental implants in humans: A comparison of three different rough surfaces. Int J Oral Maxillofac Implants 2010;25:336–344.
- Haro Adánez M, Nisihara H, Att W. A systematic review and metaanalysis on the clinical outcome of zirconia implant-restoration complex. J Prosthodont Res 2018;62:397–406.
- Wilson TG Jr. The positive relationship between excess cement and periimplant disease: A prospective clinical endoscopic study. J Periodontol 2009;80:1388–1392.
- Sailer I, Mühlemann S, Zwahlen M, Hämmerle CH, Schneider D. Cemented and screw-retained implant reconstructions: A systematic review of the survival and complication rates. Clin Oral Implants Res 2012;23(suppl 6):163–201.
- Wittneben JG, Millen C, Brägger U. Clinical performance of screw- versus cement-retained fixed implant-supported reconstructions–A systematic review. Int J Oral Maxillofac Implants 2014;29(suppl):84–98.
- Podhorsky A, Rehmann P, Wöstmann B. Tooth preparation for full-coverage restorations-a literature review. Clin Oral Investig 2015;19:959–968.
- Tiu J, Al-Amleh B, Waddell JN, Duncan WJ. Reporting numeric values of complete crowns. Part 2: Retention and resistance theories. J Prosthet Dent 2015;114:75–80.
- Reitemeier B, Schwenzer N, Ehrenfeld M, Biffar R. Einführung in die Zahnmedizin. Stuttgart: Georg Thieme Verlag KG, 2006.
- Strub JR, Türp JC, Witkowski S, Hürzeler MB, Kern M. Curriculum Prothetik. Berlin: Quintessenz Verlags, 2005.
- Hoopes W, Cushen S, DuVall N, Wajdowicz M, Brewster J, Roberts H. Failure load effect of molar axial wall height with CAD/CAM ceramic crowns with moderate occlusal convergence. J Esthet Restor Dent 2018;30:249–253.
- Simon JF, Darnell LA. Considerations for proper selection of dental cements. Compend Contin Educ Dent 2012;33:28–30, 32, 34–35; quiz 36, 38.
- Kuhn K, Rudolph H, Ahlskog L, von Koenigsmarck V, Luthardt RG. Tensile bond strength of different luting agents at the cement-restorationinterface. J Dent Res 2012;91.
- Kuhn K, Rudolph H, Ahlskog L, von Koenigsmarck V, Luthardt RG. Shortterm versus long-term shear bond strength of different luting agents. J Dent Res 2013;92 (spec. issue A):2919 [abstract].
- Abbo B, Razzoog ME, Vivas J, Sierraalta M. Resistance to dislodgement of zirconia copings cemented onto titanium abutments of different heights. J Prosthet Dent 2008;99:25–29.
- Saleh Saber F, Abolfazli N, Nuroloyuni S, et al. Effect of abutment height on retention of single cement-retained, wide- and narrow-platform implant-supported restorations. J Dent Res Dent Clin Dent Prospects 2012;6:98–102.
- Silva CEP, Soares S, Machado CM, et al. Effect of CAD/CAM abutment height and cement type on the retention of zirconia crowns. Implant Dent 2018;27:582–587.
- Tiu J, Al-Amleh B, Waddell JN, Duncan WJ. Reporting numeric values of complete crowns. Part 1: Clinical preparation parameters. J Prosthet Dent 2015;114:67–74.
- Al Hamad KQ, Al Rashdan BA, Abu-Sitta EH. The effects of height and surface roughness of abutments and the type of cement on bond strength of cement-retained implant restorations. Clin Oral Implants Res 2011;22:638–644.
- Cano-Batalla J, Soliva-Garriga J, Campillo-Funollet M, Munoz-Viveros CA, Giner-Tarrida L. Influence of abutment height and surface roughness on in vitro retention of three luting agents. Int J Oral Maxillofac Implants 2012;27:36–41.
- Kim SM, Yoon JY, Lee MH, Oh NS. The effect of resin cements and primer on retentive force of zirconia copings bonded to zirconia abutments with insufficient retention. J Adv Prosthodont 2013;5:198–203.
- Mehl C, Harder S, Shahriari A, Steiner M, Kern M. Influence of abutment height and thermocycling on retrievability of cemented implant-supported crowns. Int J Oral Maxillofac Implants 2012;27:1106–1115.

- Menon NS, Kumar GP, Jnanadev KR, Satish Babu CL, Shetty S. Assessment and comparison of retention of zirconia copings luted with different cements onto zirconia and titanium abutments: An in vitro study. J Indian Prosthodont Soc 2016;16:136–141.
- Rödiger M, Kloß J, Gersdorff N, Bürgers R, Rinke S. Removal forces of adhesively and self-adhesively luted implant-supported zirconia copings depend on abutment geometry. J Mech Behav Biomed Mater 2018;87:119–123.
- Rohr N, Brunner S, Märtin S, Fischer J. Influence of cement type and ceramic primer on retention of polymer-infiltrated ceramic crowns to a one-piece zirconia implant. J Prosthet Dent 2018;119:138–145.
- Rues S, Fugina M, Rammelsberg P, Kappel S. Cemented single crown retention on dental implants: An in vitro study. Int J Prosthodont 2017;30:133–135.
- Sadig WM, Al Harbi MW. Effects of surface conditioning on the retentiveness of titanium crowns over short implant abutments. Implant Dent 2007;16:387–396.
- Safari S, Hosseini Ghavam F, Amini P, Yaghmaei K. Effects of abutment diameter, luting agent type, and re-cementation on the retention of implant-supported CAD/CAM metal copings over short abutments. J Adv Prosthodont 2018;10:1–7.
- Kuhn K, Rudolph H, Ahlskog L, von Koenigsmarck V, Luthardt RG. Short-term versus long-term tensile bond strength to zirconia and nonprecious-alloy. J Dent Res 2014;93(spec. issue A): 1339 [abstract].
- von Koenigsmarck V, Rudolph H, Kuhn K, Ahlskog L, Luthardt RG. Influence of different cementation conditions on shear bond strength. J Dent Res 2015;94(spec. issue A):0799 [abstract].
- Blatz MB, Alvarez M, Sawyer K, Brindis M. How to bond zirconia: The APC concept. Compend Contin Educ Dent 2016;37:611–617; quiz 618.
- Gargari M, Gloria F, Napoli E, Pujia AM. Zirconia: Cementation of prosthetic restorations. Literature review. Oral Implantol (Rome) 2010;3:25–29.
- Stawarczyk B, Keul C, Eichberger M, Figge D, Edelhoff D, Lümkemann N. Three generations of zirconia: From veneered to monolithic. Part II. Quintessence Int 2017;48:441–450.
- Stawarczyk B, Keul C, Eichberger M, Figge D, Edelhoff D, Lümkemann N. Three generations of zirconia: From veneered to monolithic. Part I. Quintessence Int 2017;48:369–380.
- 44. Attia A, Kern M. Long-term resin bonding to zirconia ceramic with a new universal primer. J Prosthet Dent 2011;106:319–327.
- 45. Kern M. Bonding to oxide ceramics-laboratory testing versus clinical outcome. Dent Mater 2015;31:8–14.
- Tanış MÇ, Akay C, Karakış D. Resin cementation of zirconia ceramics with different bonding agents. Biotechnol Biotechnol Equip 2015;29:363–367.
- Oilo G, Jørgensen KD. The influence of surface roughness on the retentive ability of two dental luting cements. J Oral Rehabil 1978;5:377–389.
- Rödiger M, Rinke S, Ehret-Kleinau F, et al. Evaluation of removal forces of implant-supported zirconia copings depending on abutment geometry, luting agent and cleaning method during re-cementation. J Adv Prosthodont 2014;6:233–240.
- Rudolph H, Salmen H, Moldan M, et al. Accuracy of intraoral and extraoral digital data acquisition for dental restorations. J Appl Oral Sci 2016;24:85–94.
- Johnson GH, Lepe X, Patterson A, Schäfer O. Simplified cementation of lithium disilicate crowns: Retention with various adhesive resin cement combinations. J Prosthet Dent 2018;119:826–832.
- Wille S, Lehmann F, Kern M. Durability of resin bonding to lithium disilicate and zirconia ceramic using a self-etching primer. J Adhes Dent 2017;19:491–496.